

# RESVERATROL

Analysis Items	Specifications	Test Method
Appearance	White or off white powder	Visual
Trans-Reservatrol	98% Min.	HPLC
Emodin	0.5% Max.	HPLC
Moisture	1.0% Max.	USP
Particle Size	≥95% through 60 Mesh	USP
Heavy Metals	10ppm Max.	USP
Arsenic	2ppm Max.	AOAC
Lead	1ppm Max.	AOAC
Cadmium	1ppm Max.	AOAC
Pesticide Residual	Meet USP561	USP
<b>Microbiology</b>		
Total Bacteria	1,000 cfu/g Max.	USP
Yeast & Mold	100 cfu/g Max.	USP
E.coli	Negative	USP
Salmonella	Negative	USP
S.Aureus	Negative	USP

